




CHOICES2GO QUICK-START GUIDE




CanAssist's *Choices2Go* is an app that helps people indicate their choices in everyday life and in therapy settings. Through visual, in-the-moment "choice boards" using images from a photo library or the iPad's own camera roll, users are presented with a number of options and are encouraged to tap on their choice, which is then highlighted for several seconds.

CHARGING


Make sure your iPad is sufficiently charged before using *Choices2Go*.



TAKING PHOTOS FOR CHOICE BOARDS

1. In the iPad's home screen, tap the **Camera** icon .
2. Tap and release the shutter icon  to take a photo. Take as many photos as necessary.
3. When you're done taking photos, press and release the **Home** button .

CREATING A NEW CHOICE BOARD

1. Tap the *Choices2Go* icon  to launch the app.
2. Tap **New Board**.
3. Tap **All Photos**.
4. Tap on the photos you want to appear on your choice board.
5. Tap **Done**. Your choice board is now ready for use.
6. Tap the image that corresponds to your choice. The selected image will be highlighted in green for about four (4) seconds.
7. To indicate a different choice, tap on the desired image.

Tip: Shake the iPad to re-arrange the images on the choice board.

THE "YES/NO" BOARD

1. Tap **Yes/No** on the top left corner of the screen any time you want to access this board.
2. Tap either **YES** or **NO**. The selected option will be highlighted in green for about four (4) seconds.
3. Tap **Cancel** to return to your choice board.